

Agilia® SP PCA Agilia® SP PCA WiFi


Patient Controlled Analgesia (PCA)
syringe infusion pump



Agilia and Vigilant are registered trademarks by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information. Pictures are not contractual.



Take control of the pain with Agilia SP PCA
to deliver the prescription at your patient's pattern




Agilia Connect is as intuitive as ever
and now comes with embedded WiFi
Latest standards enable cyber secure communication

Agilia SP PCA

The Agilia pump for Patient Controlled Analgesia (PCA)

Agilia SP PCA proposes 4 PCA modes to accommodate every protocol:



Easy to use


- A detailed history available to follow your patient's events
- Intuitive workflow through a simple Menu
- Easy loading dose and clinician bolus set up
- Monitoring at a glance



Safe

- Full syringe protection with a key lock to complete a safe self-administration
- Keypad protection by password to avoid overdose of medication
- Can be used in general infusion mode with a secured hardware and software protection of the pump
- Possibility to lock the pump to the pole

Vigilant Master Med



DERS Dose Error Reduction Software

- Web application for multi-point access
- Dataset capability up to 10,000 drugs
- Configurable infusion modes



Ergonomic handset

- Comfortable shape and size for your patient
- Stroke length and soft button to get the right feeling
- Light helps patient locate the handset
- Light blinks to acknowledge PCA bolus request
- Smart compartment to store the handset